


PRODUCT DATA SHEET

Revision 1	Petri Dish	
Validfor Item-No.:	6010000/6011000/6012000/6015000/ 6015010/6015020/6015030/6013000	

1.	Description / Specification	
1.1	Description	Petri Dish, sterile 6010000: 3.5 cm 6011000: 6.0 cm 6012000: 7.0 cm 6015000: 9.0cm 6015010: 9.0cm 2 compartments 6015020: 9.0cm 3 compartments 6015030: 9.0cm 4 compartments 6013000: 15.0 cm
1.2	Dimensions	see Customer Drawing
1.3	Volume	6010000: 3ml 6011000: 6~7ml 6012000: 7~8ml 6015000: 10~40ml 6015010: 25ml in compartment 6015020: 18ml in compartment 6015030: 12ml in compartment 6013000: 25~27ml
1.4	Material / Resin	<u>Dish and Lid</u> : PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	<u>Dish and Lid</u> : clear
1.6	Sterilization	sterile
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for bacteriology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	6010000/6011000/6012000	6015000/6015010/6015020/6015030/6013000
3.1	Pieces / Bag	10	10
3.2	Pieces / Box	1000	500
3.3	Lot-No.	DDD YY XXX (day, year, product No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

